

Eutech Expert pH tester



Eutech Expert pH tester accurately measures pH. With built in ATC, an ideal solution for many applications.

Specifications

pH

Range	0.00 to 14.00
Resolution	0.1pH
Relative Accuracy	±0.1pH
Auto Calibration Points	1 - 3

Temperature

Range	0 to 50.0 °C; 32.0 to 122.0 °F
Resolution	0.1 °C / 0.1 °F
Relative Accuracy	±0.5 °C / ±0.9 °F + 1 LSD

General

Power	Includes batteries (A76 or LR44 equivalent)
Auto Off	8.5 mins after last key press
Operating Temperature	5 to 45 °C
Ingress Protection	IP67 waterproof
Dimensions (L x W x H)	170 x 45 x 30