

SERVICE & MAINTENANCE

- fume cupboard maintenance & safety checks -

FUMECARE SERVICE

To ensure the safety of fume extraction equipment, it is important that Fume Cupboard servicing is undertaken on a regular basis. This will minimise the chances of equipment failing and will also ensure that the equipment runs efficiently and, with minimal downtime and ensuring the staff using the cabinets are protected.

ABOUT US

Fumecare Limited is a family owned UK based business with access to a network of service engineers offering national coverage. All our service engineers are fully qualified and have many years of field service experience. Fumecare Limited are ISO9001-2008 Accredited and Safecontractor and CHAS Accredited in Health & Safety Standards. Ref: DZ2462



FORENSIC SERVICING & DECONTAMINATION

We specialise in the servicing of all types of Forensic Fume extraction/filtration equipment. This includes the decontamination of equipment. Emergency 24 hour call out is available on this service. We can also de-commission and re-locate equipment.

REPLACEMENT CARBON AND HEPA FILTERS

We can supply and replace HEPA and Carbon filters for most makes of fume cabinets, forensic cabinets and safety cabinets at competitive prices.

FUME CUPBOARDS

We offer a full range of standard Filtration Fume Cupboards from 600mm wide up to 1700mm wide. We also have Downflow Work Stations in three sizes 600mm wide, 900mm wide and 1200mm wide. We also supply a range of Forensic Equipment such as Forensic Drying Cabinets, Fingerprint Benches and Superglue Chambers. We also offer a full custom design service.

FUME CARE