

# Haier Biomedical

## Laboratory Refrigerator

Adjustable Temperature Set Point

Precise Temperature Control

Designed specifically for labs



\* Temperature control is precise and stable



HLR-198F/310F

- **Wide Range Temperature Set Point:**  
3-16°C Adjustable set point with one key – flexible and adaptable to meet your needs.
- **Precise Temperature Control:**  
Internal temperature variation is within 3°C to professionally keep temperature uniformity.
- **Energy Saving and Environmentally Friendly:**  
Hydro carbon cooling system and energy saving fan motor provide environment friendliness and energy saving.
- **Multiple Alarms:**  
Include high and low temperature alarms, door open alarm, etc. The refrigerator can be remotely monitored and will alert users rapidly.



Scan 2D barcode with your mobile device to learn more information about Haier Biomedical

### Haier Biomedical UK Ltd.

Address: Mytogen House, 11 Browning Road, Heathfield, TN21 8DB, United Kingdom

Tel: (01435) 517 003

Fax: (01825) 700 471

Email: [sales@haierbiomedical.co.uk](mailto:sales@haierbiomedical.co.uk)

[www.haierbiomedical.co.uk](http://www.haierbiomedical.co.uk)

**Haier Biomedical**  
Intelligent Protection of Life Science

**Main Uses** ➤ To store chemical and experiment reagents and consumables for laboratories, research facilities, etc.

### Product Advantages



#### 3-16°C wide temperature range with precise control

3-16°C wide range temperature adjustment. Intelligent temperature control system keeps internal temperature variation within 3°C for great temperature uniformity.



#### One Key Search for Convenient Checking

Historical temperature data can be saved and checked by pressing only one key for convenient laboratory inspection.



#### Dual Lock Design with Multiple Alarm System

The door comes with a mortise lock. A padlock can also be installed. The system has high and low temperature, open door, and sensor error alerting functions for better safety.

### Ergonomic Design



#### HC Cooling Saves Half Energy

The hydro carbon cooling system has twice the cooling efficiency and saves half the energy.



#### Low Noise

Ultra quiet fan motor with noise cancelling technology makes the entire refrigerator's noise lower than 40 decibels for a quiet and comfortable lab environment.



#### Glass Shelves

Tempered glass shelves can store over 40kg of various sized reagent bottles.



#### Automatic Door Closing Function

Provides small angle automatic door closing function.



Cabinet Structure



Glass Shelves



Fan/Interior Light



Display Board

### Product Specifications

Model	Voltage	Power	Internal Temperature Range	Exterior Dimensions (W*D*H)	Interior Dimensions (W*D*H)	Capacity	Shelves	Net Weight
HLR-198F	220V/50Hz	195W	3°C~16°C	535*570*1530	455*445*900	198L	4	54kg
HLR-310F	220V/50Hz	270W	3°C~16°C	605*598*1840	560*500*1290	310L	4	77kg