



Laboratory & Mortuary Equipment

CO₂ Incubator Accessories




Model PRV



Model PNEU



Model R06



A typical CO₂ installation

Model	Description	Price
R06	In-line CO ₂ reducing valve with pressure gauge (2 to 30 p.s.i.)	£49.00
PRV	Single-stage CO ₂ cylinder pressure reducing valve with gauges	£77.00
PRV2	Two-stage CO ₂ cylinder pressure reducing valve with gauges	£139.00
PRVN	N ₂ cylinder pressure reducing valve with gauges	£139.00
PNEU	Automatic change over unit for two cylinders (CO ₂ or N ₂)	£742.00