



AUTOCLAVES

MDB18C

Table top autoclaves for the modern laboratory, the pharmaceutical and food industry as well as for special applications. Innovations for safe, accurate, reproducible and validation for sterilisation in research and industry within the dental, medical, cosmetic, laboratory, veterinary and tattoo fields.

Our autoclaves are produced under the most rigorous system of quality control certifications in accordance to ISO13485 and TUV CE0197.

All models are fully automatic with LCD display, independent steam generator, temperature and pressure monitoring systems. Every unit comes equipped with a unique high precision control system and automatic alarm for water level, door locking and cycle errors.

Product Code	MDB8C	MDB18C	MDB23C
Volume(Litres)	8	18	23
Chamber Size(mm)	170x320	247x350	247x450
Power (W)	1500	2200	2200
Printer (external)	Optional		
Double Door locking	Yes		
No.of trays (max)	2	4	
Dimensions (mm) W X D X H	370x565x345	480x700x455	480x700x455
Net Weight (Kg)	34.5	51	58

OPTIONAL:

- Printer
- USB port and USB stick
- Double door locking device
- Independent Drying programme



MDB23C



MDB8C

PROGRAM	TEMP (°C)	PRESSURE (Mpa)	STERILIZATION TIME (Min)	TOTAL TIME		
				B8C	B16C & B23T	B18C & B23C
SOLID	134	210	4	15-25	14-30	20-45
	121	110	20	25-40	30-45	35-60
LIQUID	134	210	10	25-50	30-55	35-60
	121	110	30	30-55	35-60	35-60
WRAPPED	134	210	4	20-45	30-45	45-65
	121	110	20	30-50	30-50	50-75
TEXTILE	134	210	4	20-45	30-45	45-65
	121	110	20	30-50	30-50	50-75
PRION	134	210	18	30-50	35-55	45-75
DRYING	-----	-----	-----	1-20	1-20	1-20
B&D TEST	134	210	3.5	22-35	22-35	22-35
HELIX TEST	134	110	3.5	22-35	22-35	22-35
VACUUM TEST	-----	-----	-----	15-20	15-20	15-20