

Data Handling Made Direct and Easy With EasyDirect™ pH Software



Simplicity of data handling

Using EasyDirect™ is so easy that there is no need to learn how to use or reading its manual. It is intuitive to use and feels familiar in its operation, just what you would expect. Additionally the application is available in 12 languages.



Plug and Play connectivity

Complex communication settings are past with EasyDirect pH, simply connect your SevenCompact™, Seven2Go™ or FiveEasyPlus to your computer, confirm that you like to add this instrument and start using the full functionality of the software.



Intuitive operation

To start immediately the software is set to display most common parameters of your instrument by default. Nevertheless you can modify within seconds the displayed parameters according to your needs. For your convenience you can save different screen layouts to call up display, export or printing.



Data Handling Made Easy

Connect your pH meter to EasyDirect™ pH software to improve your data management. Within a few clicks the data measured on your pH meter could be seamlessly transferred to your computer. The application will store the data unmodified and allow for further printing or analysis. Enjoy the possibility to export the data as you see them to excel or any other csv application. Directly print the data to any printer connected or installed on your computer. Relax in the knowledge that your measurement results are stored and easily accessible with EasyDirect™ for further analysis.

Download your 45 days free trial version on:

► www.mt.com/EasyDirectpH

Seamless communication

Need the data for statistical analysis, or do you need a physical printout for archiving? Within few seconds EasyDirect™ pH can export the data unmodified to Excel or any other column separated application. Also printing to any printer connected to your computer is simple. The software will print the selected data as you see them and fit the column width to one page.



Connectivity

Connect SevenCompact™, SevenCompact™ Duo, Seven2Go™ Pro, FiveEasyPlus or SevenExcellence™ to EasyDirect™ pH Software either by USB or RS232 connection. Up to three instruments could be connected simultaneously to EasyDirect™ pH and up to three instruments data can be stored.



Mettler-Toledo GmbH, Analytical

Im Langacher 44
8606 Greifensee, Switzerland
Tel. +41 22 567 53 22
Fax +41 22 567 53 23

Subject to technical changes
© 06/2016 Mettler-Toledo GmbH, 30353939
Marketing pH Lab / Global MarCom 1957 MA

www.mt.com/EasyDirectpH

For more information