

## PM-6 ADVANCED EDUCATIONAL MICROSCOPE WITH LED ILLUMINATION



The PM-6 microscope is an economically priced instrument but with a performance for exacting microscopy. The focusing Abbe condenser provides suitable illumination for use with the 100x oil objective. The built in long life LED illumination provides bright even illumination.

The objectives are DIN corrected. Ergonomically positioned controls, including drop down mechanical stage controls provide comfortable viewing.

A version is available with a binocular head with phototube enabling a camera to be readily attached.

Educational Microscopes		
Product:	Description:	Price:
PM-6	<ul style="list-style-type: none"> <li>• Binocular head with widefield 10x/18 eyepieces</li> <li>• Quadruple nosepiece</li> <li>• Objectives 4x, 10x, 40x, 100x oil</li> <li>• Built-in drop control mechanical stage</li> <li>• Focusing Abbe condenser with iris</li> <li>• Coaxial coarse &amp; fine focusing</li> <li>• Built-in 3W LED illumination</li> </ul>	£ 670.00
PM-6TR	As PM-6 except with a trinocular head, enabling a camera to be fitted	£ 740.00

## Additional Components

Extra Objectives		
Product:	Description:	Price:
PM-O20P	Objective 20x	£ 50.00
PM-O60P	Objective 60x	£ 100.00

Measuring		
Product:	Description:	Price:
PM-E10FG	Wide-field focusing eyepiece 10x/20 with scale graticule 10:100	£ 60.00
PM-STG	Calibrating stage micrometer 1mm in 100 parts	£ 40.00

Camera Adapters		
Product:	Description:	Price:
PM-6CC1	C-mount adapter 0.5x	£ 100.00
PM-6CC2	C-mount adapter 0.3x	£ 100.00
LM-PT	Phototube	£ 60.00

### Listed prices exclude delivery & VAT

TEC Microscopes LTD, The Vineyard, Axbridge, Somerset, BS26 2AN  
 Tel: 01934 733679 | Fax: 01934 733660  
 Email: [enquiries@tecmicroscopes.co.uk](mailto:enquiries@tecmicroscopes.co.uk)