



More information on the website  
radwag.com/en/info,w1,DQV

## WTC 3000 Precision Balance



## Functions



Plus/Minus Control



Percent Weighing



Parts counting



Peak hold



GLP Procedures



ALIBI Memory

## Datasheet

| Metrological parameters |               |
|-------------------------|---------------|
| Maximum capacity [Max]  | 3100 g        |
| Readability [d]         | 0,1 g         |
| Tare range              | -3100 g       |
| Repeatability           | 0,1 g         |
| Linearity               | ±0,3 g        |
| Stabilization time      | 2 s           |
| Adjustment              | external      |
| Physical parameters     |               |
| Display                 | LCD (backlit) |

| Physical parameters         |                                       |
|-----------------------------|---------------------------------------|
| Weighing pan dimensions     | 128×128 mm                            |
| Packaging dimensions        | 330×220×140 mm                        |
| Net weight                  | 1,3 kg                                |
| Gross weight                | 2 kg                                  |
| Protection class            | IP 43                                 |
| Communication interface     |                                       |
| Communication interface     | RS232, USB-A, USB-B, Wi-Fi (option)   |
| Electrical parameters       |                                       |
| Power supply                | 100 ÷ 240 V AC 50 / 60 Hz and battery |
| Power consumption max.      | 6 W                                   |
| Operation time on batteries | 15 h (average time)                   |
| Environmental conditions    |                                       |
| Operating temperature       | +15 ÷ +30 °C                          |

\* Wi-Fi® is a registered trademark of Wi-Fi® Alliance.



## Accessories

Receipt Printer  
Cigarette lighter receptacle power supply cables



RS 232, RS 485 cables  
RS 232 cables (scale - EPSON printer)

## Software

RAD-KEY  
LabVIEW Driver  
R-LAB

RW Tool  
R Panel  
Scales Editor 2.0

## Device dimensions



WTC: d = 0.01 g, d = 0.1 g