



HOTPLATES

**MODEL II SERIES—GENERAL PURPOSE
P.T.F.E. COATED, ALUMINIUM OR
CAST IRON TOP OPTIONS.**



GENERAL DESCRIPTION.

This range of hotplates has been specifically designed for use in the general laboratory. The robust top casting, which is heated by an incoloy sheathed rod element, has angled sides to discard accidental spillage away from the base plinth which houses the temperature control and other electrical components. Two red baffle plates are positioned between the top casting and the base plinth to give extra protection from spillage and heat transfer. Temperature control is either by simmerstat or thermostat

PRICES: are exworks. Packing & Delivery charges are for Buyer's account & subject to alteration without notice.

FAX: 01279-438095.

TEL: 01279-422456.

PLEASE SEE PAGE No. 4. FOR TYPES, SIZES & PRICES.

