



HOTPLATES



CERAMIC TOP



GENERAL FEATURES:

- * Easy Clean, No Maintenance Ceramic Top
- * Acid & Alkali Resistant
- * Rapid Heatup
- * Maximum Temperature 450 C
- * Hot Zone Indication
- * Constructed from Non Rusting Materials

S & J. JUNIPER & CO.

7, POTTER STREET, HARLOW. ESSEX, U.K.

Telephone:- 01279 422456.

Fax:- 01279 438095.

Email: info@sjjuniper.com Web: www.sjjuniper.COM

