



HOTPLATE



MODEL H116

PRICE: GBP 739.25.

Robust one piece cast iron construction, Can be free standing on a bench or recessed into a fume cupboard to give a flush heated area.

- **General Description.**
- The unit is constructed around a one piece cast iron casting having a machined flat top surface . Incaloy sheathed rod elements are clamped to the underside of the top surface with a downward heat radiation being limited by ceramic fibre thermal blanket which is held in position by a zinc coated mild steel tray.
- Power cables and thermocouple leads are housed in a flexible housing which is connected to the **Digital** controller unit.

Temperature Control.(Now Digital.)

This is achieved by the use of a “type K” thermocouple, mounted in direct contact with the underside of the heated surface and linked to a **Digital** temperature control unit incorporating “MAINS” and “LOAD” indicators.

MAXIMUM TEMPERATURE: APPROX. 410 Deg C.

SUPPLY VOLTS: 230/240 V.A.C.

TOP PLATE DIMENSION: 457 mm.SQUARE.

LOADING: 3.9 Kw.

S & J.JUNIPER & CO.
7, POTTER STREET, HARLOW, ESSEX, CM17 9AD. U.K.
TEL: 01279-422456 FAX: 01279-438095.
WEB: www.sjjuniper.com EMAIL: info@sjjuniper.com

