



LABORATORY HOTPLATES ROUND SERIES ALUMINIUM TOPS.



H111. PRICE: GBP 180.54.

THIS UNIT IS NOT FITTED
WITH A MEANS OF
TEMPERATURE CONTROL



TOP.

Cast Aluminium 150 mm Dia.
With ground top surface.

BASE.

Cast Aluminium with three
Integral cast feet to stand on
bench and a screw "Vee" clamp
For use on a retort stand.

PERFORMANCE.

Temperature rise of
HotplateTop
From ambient to 415 deg C.
In approx. 15 min.
Maximum temperature
of approx 450 deg.C.

TOP.

Cast Aluminium 200 mm Dia.
With ground top surface.

BASE.

Cast Aluminium with three
integral cast feet to stand on
bench and a screw "Vee" clamp
for use on a retort stand.

BOTH 230/250 V.A.C.

H112S. PRICE: GBP 229.64.

H112T. PRICE: GBP 239.49.

P.T.F.E. COATED.

H112SP. PRICE: GBP 262.96.

H112TP. PRICE: GBP 272.50.



CAT.No. TEMP.

**LOADING APPROX.TEMP.RANGE.
CONTROL**

H112/S SIMMERSTAT. 400 WATTS 50 - 420 Deg. C.

S & J.JUNIPER & CO, 7, POTTER STREET, HARLOW, ESSEX. CM17 9AD.

TEL: 01279-422456.
WEB: www.sjjuniper.com



FAX: 01279-438095.
EMAIL: info@sjjuniper.com