

## SF RANGE OF -20°C FREEZERS

| MODEL                      | SW425 SFW        | SW625 SFW        | SW660 SFS        | SW1270 SFS        | SW1400 SFS        |
|----------------------------|------------------|------------------|------------------|-------------------|-------------------|
| EXT DIMS (mm)<br>H x W x D | 2000 x 600 x 735 | 2000 x 815 x 735 | 2150 x 815 x 760 | 2000 x 1390 x 875 | 2150 x 1390 x 875 |
| CAPACITY (litres)          | 425              | 625              | 660              | 1270              | 1400              |
| DOORS (shelves)            | 1 (4)            | 1 (4)            | 1 (5)            | 2 (8)             | 2 (10)            |
| WEIGHT (kgs)               | 115              | 125              | 130              | 250               | 260               |

- Adjustable temperature range -15°C to -25°C
- Fan circulation for even temperature distribution
- Adjustable wire shelves
- Lockable door(s)
- SFW models have white exterior
- SFS models have stainless steel exterior
- Electronic control
- Digital temperature display
- Audio/visual temperature alarm
- Automatic defrost
- Option: Access port for monitoring probe
- Option: Stainless steel exterior (for SFW models)

These models are not internally sparkfree

SW625 SFW



SW660 SFS



SW1270 SFS

