

VFE RANGE OF -20°C FREEZERS

MODEL	SW200 VFE	SW400 VFE	SW550 VFE	SW600 VFE
EXT DIMS (mm) H x W x D	850 x 600 x 600	1850 x 600 x 600	1720 x 780 x 720	1890 x 780 x 735
CAPACITY (litres)	140	400	550	600
SHELVES	2 Fixed	7 Fixed	7 Fixed	7 Fixed
WEIGHT (kgs)	45	75	95	105

- Adjustable temperature range -18°C to -25°C
- Electronic control
- Digital temperature display
- Internally sparkfree
- Visual temperature alarm
- Static cooling (fixed freezing shelves)
- Lockable door
- Manual defrost
- White exterior
- Stainless steel exterior available
- Options: Audible alarm, access port for monitoring probe

SW200 VFE
(white)



SW400 VFE
(stainless steel)

