

Vented Storage Cabinets

A powerful fan at the top of the cabinet draws air at 90 cubic metres/hour throughout the enclosed cabinet passing it through a filter which removes noxious pollutants and odours. The purified air is then returned to the room. Cabinets are manufactured from a chemical resistant alloy and stainless steel, coated with white acid-resistant paint. Available with shelves or sliding drawers.

- Also purifies room odours
- Consumes just 47 Watts
- Spill tray located at base of cabinet
- Sparkless free motor will not ignite fumes
- Levelling jacks compensate for uneven floors

With Shelves

Product No	11501/01
Model	AVP 804
No of shelves	One fixed Four adjustable
Door	Lockable plus magnetic fastening. Made of clear self- extinguishing PMMA. Opens 180°
External Dimensions	2150mm high 490mm deep 800mm wide
Storage	100 x 1Lt bottles

With Sliding Drawers

Product No	11501/02
Model	AVVE 804
No of drawers	Three. Upper drawer is lockable
Sliding mechanism	Ball bearing slides
External Dimensions	1640mm high 600mm deep 800mm wide
Storage	100 x 1Lt bottles & 400 x 250ml bottles
Max storage weight	50kg per drawer 150kg total

