

Disposal Boxes

Broken Glass

Sturdy, disposable boxes of corrugated cardboard. A 50 micron thick inner polyethylene bag retains the glass together with associated liquids. The close fitting lid has a flap through which waste material is safely deposited. When full, the safety cap is pulled into place making the entire unit ready for disposal. Two sizes are available for benchtop and floor use. Supplied flat to save storage space, they take only seconds to assemble. Supplied in packs of 6



Product No	Size (H x W x D)-inches
------------	-------------------------

11353/01	27 x 12 x 12
----------	--------------

11353/02	10 x 8 x 8
----------	------------

Biohazards

Similar to above but clearly marked "Biohazard Incinerator Carton", with a 40 micron HMHD biohazard bag. Supplied in packs of 6

Product No	Size (H x W x D)-inches
------------	-------------------------

11354/01	27 x 12 x 12
----------	--------------

11354/02	10 x 8 x 8
----------	------------

Safety Waste System

Solvent, Chemical & Biological Waste

Each all-plastic system consists of an H.D.P.E funnel with cover and waste container to reduce volatile emissions. Suitable for temporary benchtop storage, a built-in vent releases excess pressure



Product No	Capacity-Lt	Funnel Dia. Top I.D.-mm
------------	-------------	-------------------------

11355/01	2	102
----------	---	-----

11355/02	4	102
----------	---	-----

11355/03	10	203
----------	----	-----

Accessories

It is recommended that for additional safety, the above systems are placed within one of the following secondary containers to contain any spills. (11355/13 is supplied with a lid).

Product No	Capacity-Lt	Material	To suit
11355/11	3.9	LDPE	11355/01
11355/12	5.5	LDPE	11355/02
11355/13	18	HDPE	11355/03

Smoke Alarm

Ionisation type smoke alarm. Battery powered with test button. (Battery should be ordered separately).

- 85dB piezo-electric alarm
- Dual ionisation chambers
- Easy to install
- LED operating light
- Low battery warning

Product No	Description
------------	-------------

11279/01	Smoke alarm
----------	-------------

10993/05	9 volt alkaline battery (PP3)
----------	-------------------------------

Exit Signs

Illuminated 8 watt fluorescent exit sign in white epoxy coated sheet steel housing with rear, top or side cable entry. Measures 195mm high, 425mm wide, 65mm deep. Sign illuminates for 3 hours in the event of a power failure. Choice of diagrams to be ordered separately

Product No	Size
------------	------

11356/01	Sign housing
----------	--------------

11356/02	Sign legend, arrow down
----------	-------------------------

11356/03	Sign legend, arrow right
----------	--------------------------

11356/04	Sign legend, arrow left
----------	-------------------------

Heat Resistant Mats

These mats provide heat protection to benches. Made of 5mm thick ceramic fibre, they are a safe alternative to asbestos

Product No	Size-mm
11357/01	150 x 150
11357/02	225 x 225
11357/03	300 x 300

