



...the technical stuff

Features

- 0.2µm Sartorius PTFE filter
- Filter self-sterilisation
- Cooling fan
- Suitable for BSL-3 labs as supplied
- 316L stainless steel chamber
- Complies with ASME, CE, DOHS, MOM and TS
- Panel mounted printer for temperature, time and F0 value.
- Digital pressure sensor
- Load sensed process timing with wandering PT100 probe
- Saturated steam monitoring
- Locking castors
- Auto drying (SRU model only)

Weights and dimensions:

Model	AMC-GR60-SFU	AMC-GR60-SRU	AMC-GR85-SFU
Capacity	60 litres		85 litres
Chamber dimensions	390mm dia x 505mm deep		390mm dia x 700mm deep
Basket dimensions	370mm dia x 190mm deep		370mm dia x 280mm deep
Basket capacity	2		
Overall dimensions	600x831x980mm (LxWxH)		644x831x980mm (LxWxH)
Weight, dry, empty	111kg		120kg

Power requirements:

Model	AMC-GR60-SFU	AMC-GR60-SRU	AMC-GR85-SFU
Rated power	2.9kW		4.6kW
Power requirement	220V±10%, 16A, 50/60Hz		220V±10%, 32A, 50/60Hz

Performance:

Model	AMC-GR60-SFU	AMC-GR60-SRU	AMC-GR85-SFU
Temperature display accuracy	0.1°C		
Working environment	5-40°C, RH ≤85%		
Chamber material	316L		
Sterilisation temperature range	105-138°C		
Sterilisation time range	1-6,000min (100 hrs)		
Melting temperature range	60-115°C		
Melting time range	1-6,000min (100 hrs)		
Warming temperature range	45-79°C		
Warming time range	1-9,999min (approx. 167 hrs)		
Drying time		1-300min (5 hrs)	
Drying temperature		80-160°C	
Cooling lock range	40-99°C for solids, discard, agar 40-80°C for liquids		
Waste water tank	3 litres		
Cooling water tank		6 litres	
Sterilisation programmes	Liquids, Liquids with warming, Rubber, Fast, Discard, Agar, +User defined	Liquids, Liquids with warming, Rubber, Fast, Discard, Agar, Wrapped instruments, Fabrics, +User defined	Liquids, Liquids with warming, Rubber, Fast, Discard, Agar, +User defined
Programme storage	60 programmes		
Exhaust speed control	Fully automatic internal discharge, choice of 6 levels		
Time delay for auto start up	0-15 days		
Pressure gauge range	-0.1 – 0.5MPa		
Safety valve set pressure	0.29MPa		
Controller	"Inspiration" controller with LCD display		
Lid opening	Easy one-hand operation		
Locking device	Patented 8-point locking		

Reagent bottle capacities:

Model	AMC-GR60-SFU/SRU		AMC-GR85-SFU	
Bottle capacity	No. of baskets	Total No. of bottles	No. of baskets	Total No. of bottles
250ml	2	36	2	38
500ml	2	24	2	24
1000ml	1	8	2	16
2000ml	1	4	2	8
5000ml	1	1	1	1
10000ml			1	1

Erlenmeyer flask capacities:

Model	AMC-GR60-SFU/SRU		AMC-GR85-SFU	
Flask capacity	No. of baskets	Total No. of flasks	No. of baskets	Total No. of flasks
250ml	2	26	2	26
500ml	1	7	2	14
1000ml	1	5	2	10
2000ml	1	2	1	2
5000ml			1	1

Accessories:

AMC-GR60-BASKET

Autoclave Basket, 190mm deep (for 60 litre model)

AMC-GR85-BASKET

Autoclave Basket, 280mm deep (for 85 litre model)

amicus

Due to our policy of continuing product improvement the above specifications are for guidance only and subject to change without notice.